

rematitan[®] CL  powered by Dentaureum

Titanium alloy (powder) acc. to EN ISO 9693/DIN EN ISO 22674, Type 4

CE 0483

With an appropriate approval* rematitan[®] CL can be used for production of metallic restorations by means of the metal laser melting process.

22

Ti

47,88

RANGE OF APPLICATION

With an appropriate approval* rematitan[®] CL can be used for production of crowns and bridges, frames for metal ceramic veneering, cast partials, primary – and secondary parts for combined restorations as well as for implant supraconstructions.

CHEMICAL COMPOSITION

Component	Mass (%)
Ti	90
Al	6
V	4
Other elements <1 %: N, C, H, Fe, O	

TECHNICAL DATA IN LINE WITH DIN EN ISO 9693 / DIN EN ISO 22674 AFTER RECOMMENDED HEAT TREATMENT


Yield Strength $R_{p0,2}$	950 MPa
Tensile Strength R_m	1005 MPa
Elongation at fracture A_5	10 %
Modulus of elasticity E	115.000 MPa
Melting range Δ	1604-1655°C
Density ρ	4,5 g/cm ³
Coefficient of thermal expansion TEC (25-500°C)	10,16 x 10 ⁻⁶ K ⁻¹
Colour	white
Metal- ceramic bond strength acc. to EN ISO 9693, 3-Pt.- bending test (min. 25 MPa acc. to EN ISO 9693)	37 MPa (Triceram, Dentaum)
Type	4
Biocompatibility, L 929-Proliferation acc. to EN ISO 10993-5, -12	No deliberation of cell toxic active substances
Corrosion resistance, static immersion test acc. to EN ISO 10271 (max. 200 µg/cm ² x 7d acc. to EN ISO 22674)	Ion release 1,41 µg/cm ² x 7d

rematitan[®]
CL


powered
by Dentaum

MICROSECTION

Test piece etched (20 x)



Test piece etched (100 x)



HEAT TREATMENT

Perform heat treatment under an argon atmosphere. Heat up within 4 hours to 820°C. Maintain temperature for 1,5 hours. Allow the components to cool down to 500°C in the oven.

MICROSTRUCTURE

Components made from the titanium alloy rematitan[®] CL display a homogeneous, pore-free structure after they are constructed by means of the metal laser melting process LaserCUSING[®].

Concept Laser GmbH
 An der Zeil 8
 D 96215 Lichtenfels

T: +49 (0)95 71. 1679 200
 F: +49 (0)95 71. 1679 299
 info@concept-laser.de

Concept Laser Inc.
 1000 Texan Trail, Ste 150
 Grapevine, TX 76051 USA

T: +1-817-328-6500
 info@conceptlaserinc.com
 www.conceptlaserinc.com

All of the specified figures are approximate figures. The figures which are provided reflect the current level of our knowledge and are dependent on process and machine parameters. The information provided on this material data sheet is therefore not binding and is not deemed to be certified.
 * The approval is branch-specific and/or application-specific and it must be, therefore, carried out by the consumer/user. Approval of materials by Concept Laser GmbH is not available.