



PRODUCT INSTRUCTION SHEET



REFERENCE Hallaton Mobile Shower Chair (141)
141-2 / 141-4 / 141D-2 / 141D-4

DESCRIPTION

Tubular metal frame with a moulded polypropylene
Seat with small drainage holes. The arms are available
Fixed or detachable and the shower chair is mobile
With 2 or 4 foot operated brake castors and a pushing
Handle.

TECHNICAL SPECIFICATION

Seat width x depth:	18" x 15"
Fixed Seat height:	19"
Outside depth x width	20" x 21"
Internal width between arms:	19"
Maximum user weight:	200kg (31.5st)
Product weight:	8kg (18lb)



SAFETY WARNING

Ensure that at least two of the castors are braked when the shower chair is stationary, particularly when a person is getting in and out of it.

CLEANING AND CARE

To clean, use a mild detergent, and warm water. Allow all equipment to dry thoroughly.

WARRANTY

Your Welland by Drive branded product is warranted to be free of defects in materials and workmanship for one year from date of purchase. This device was built to exacting standards and carefully inspected prior to shipment. This warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the customer. In the event of a defect covered by this warranty, we will at our option supply parts or replace the device. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, castors and grips, which are subject to normal wear and need periodic replacement. If you have any queries or require full warranty conditions, please contact your Drive stockist. This does not affect your statutory rights.